


www.genscopharma.com

info@genscopharma.com

@



# Flowflex SARS-CoV-2 Antigen Rapid Test

Qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in anterior nasal swab specimens directly from individuals within 7 days of symptom onset or without symptoms or other epidemiological reasons to suspect COVID-19 infection. This test is authorized for non-prescription home use with self-collected anterior nasal swab specimens directly from individuals aged 14 years and older or with adult-collected anterior nasal samples directly from individuals aged 2 years or older. **The Flowflex COVID-19 Antigen Home Test does not require serial testing.** 


- Anterior nasal swab specimens
- Results in 15 minutes
- 12 Months shelf life
- Store between 36 to 86° F

- · Sample self-collection ages 14 and older
- Sample by an adult in children ages 2 to 13
- Excellent performance when compared to an FDA authorized molecular SARS-CoV-2 test.

### **Clinical Performance**

The Flowflex COVID-19 Antigen Home Test was compared to an FDA authorized molecular SARS-CoV-2 test. The Flowflex COVID-19 Antigen Home Test correctly identified 93% of positive specimens and 100% of negative specimens.

### Materials Provided



### Ordering Information

Product Name	Catalog No.	Format	Specimen	Package
Flowflex COVID - 19 antigen Home Test	L031-118B5	Cassette	Nasal Swabs	1 Test/Kit
Flowflex COVID - 19 antigen Home Test	L031-125P5	Cassette	Nasal Swabs	25 Tests/Kit

This product has not been FDA cleared or approved but has been authorized by FDA under an EUA.

• This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.


• The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of IVDs for detection and/or diagnosis of

COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated, or authorization is revoked sooner.


· For more information on EUAs please visit: https://www.fda.gov/emergency-preparednessand-response/mcm-legal-regulatory-and-policy-framework/emergency-useauthorization

• For the most up to date information on COVID-19, please visit: www.cdc.gov/COVID19













# **Manufacturer SKU** UPC



Tests Kits/Sale Pack Sale Pack Dims (LxWxH) Sale Pack Weight (lbs) Sale Pack Cube





Cartons/Pallet Sale Packs/Pallet Test Kits/Pallet TI-HI Pallet Dims (LxWxH) Pallet Weight (lbs)

Sale Packs/Carton

Carton Dims (LxWxH)

Carton Weight (lbs

Test Kits/Carton



Pallets/Truckload Sale Packs/Truckload Test Kits/Truckload Truckload Weight (lbs)

# 1 Test/Kit ACON-L031-118B5 68260766026

1 4.72"x2.55"x1.06" 0.07 0.0073

300 300 22.24"x15.35"x14.37" 28.00

> 25 7,500 7,500 5x5 40"x48"x73" 740 24 180,000 180,000

> > 17,760

@

### 25 Tests/Kit L031-118B5 68260766026

25 8.27"x4.92"x3.66" 0.85 0.0861

24 600 17.51"x17.32"x16.14" 20.29

> 24 576 14,400 6x4 40"x48"x " 487

# 24 13,824 345,600 11,687



